



UNIVERSITY OF
LATVIA

Visible-NIR Multispectral Imaging of Whole-Body Human Skin

Janis Spigulis, Edgars Kviesis-Kipge, Uldis Rubins, Ilze Irbe, Marta Skrastina, Jekaterina Tihomirova

Biophotonics Laboratory, Institute of Atomic Physics and Spectroscopy
Faculty of Science and Technology, University of Latvia

NBC-20, Warsaw (PL), 17/06/2025

Whole-body skin imaging



Canfield Vectra WB360 3D




← Multi-camera design

Single camera design →



This study:
broadband white illumination
replaced by multi-laser line
illumination → extracting sets
of **narrowband spectral
skin images** in the visible
and near-infrared ranges for
improved diagnostics





Fotofinder ATBM master 4th Gen



The main idea



Combined spectral lines:
450nm, 520nm, 638nm (VIS),
+ 850nm, 940nm (NIR)
→ subsequent extraction of
skin spectral images at these
wavelengths





Skin penetration/imaging
depth increases with
wavelength up to ~1100 nm



White broadband illumination:
halogen lamp, LED or Xe lamp

Extraction of 3 spectral images at triple laser line illumination

- Three spectral images, each related to one specific wavelength, can be extracted from a single-snapshot image data. Conditions: the relative sensitivities of photo-detector at all three illumination wavelengths are known and the photo-response is linear (Proc.SPIE, 2014)
- Illumination combinations exploited: (i) 450/520/638nm; (ii) 450/520/850nm; (iii) only 940nm



Equipment: prototype system



61-megapixel color camera (Sony a7R IVA), NIR filter removed



Side-emitting optical fiber spirals: 60 m (RGB laser) and 30m (850nm laser)

3 x 1 W RGB laser with SMA output, 450/520/638 nm + 1W 850nm + 940nm LED set

Camera & illuminator up-down movements (0.5m - 1.5 m)

Placed in a light-shielding tent or in a dark room (Latvian Oncology Center patient measurements)

Clinical study in LOC: 42 patients, 63 malformations




Previous research:


J.Spigulis, U.Rubins, E.Kviesis-Kipge, I.Saknite, I.Oshina, E.Vasilisina, **“Triple spectral line imaging of whole-body human skin: equipment, image processing, and clinical data”**, *Sensors* **24**, 7348 (2024).

DOI:10.3390/s24227348.

Distribution of lesion location on the body, %



The image processing tool



Examples of clinical spectral images: MM on the back skin


Color/RGB image

450 nm

520 nm

638 nm

850 nm



Examples of clinical spectral images: BCC (basalioma)

Color/RGB image



450 nm

520 nm

638 nm

850 nm

940 nm




Provisional results: contrasts in malformation's spectral images




$$\text{Contrast} = I_{\text{median_skin}} / I_{\text{median_malf}}$$

$I_{\text{median_skin}}$ - median pixel value of surrounding skin

$I_{\text{median_malf}}$ - median pixel value of malformation



Provisional results: standard deviations in histograms of spectral images



Summary

- A prototype system for multispectral imaging of full-body skin at five VIS-NIR wavelengths (450/520/638/850/940 nm) was assembled and clinically tested
- 63 skin malformations were examined in Latvian Oncology Center, including 50 malignant tumors - 4 melanomas and 46 basal cell carcinomas
- General trend: most skin malignancies are still visible in NIR spectral images, while nevi and other benign malformations fade/disappear. Could be helpful for fast early screening of skin melanoma.
- Potential diagnostic criteria for MM, BCC and nevi: contrasts and histogram standard deviations yielded from NIR spectral images, or their ratio to those at visible wavelengths (e.g. 520nm/850nm)
- More patient data should be collected, our clinical study continues

Acknowledgments

This research is funded by the Recovery and Resilience Facility project "Internal and External Consolidation of the University of Latvia" (No.5.2.1.1.i.0/2/24/I/CFLA/007), grant No. [LU-BA-PA-2024/1-0006](#).



**Funded by
the European Union**
NextGenerationEU



Thank You!

- [*janispi@latnet.lv*](mailto:janispi@latnet.lv)
- [*http://home.lu.lv/~spigulis*](http://home.lu.lv/~spigulis)

